



~ SIREN Elite Training Highlight ~

HOW TO RECORD A 12-LEAD IN SIREN ELITE

In the old SIREN, you could record a 12-Lead in the procedures list, or in the “Add EKG” section within the Vitals/Treatment tab. In the new SIREN Elite, the 12-Leads have been taken out of the procedure list, so you will record all your 12-leads in the ECG grids on the Treatment Panel, located within the Provider Action section.



There are two ECG grids on this treatment tab:



This one is for services who have their cardiac monitors linked to their run form

If your cardiac monitor is linked with your run form, all your 12-lead data, as well as defibrillation and vitals information collected by the monitor, will appear in this “AED/Cardiac Monitor Imports” grid.

AED / Cardiac Monitor Imports

+ Add Another **✓ OK** **✕ Cancel**

Date/Time of Event (Captured by Medical Device):

Medical Device Event Type:

Medical Device ECG Lead:

Medical Device ECG Interpretation:

ECG Cardiac Rhythm:

EKG Cause for Change:

Total Number of Shocks Delivered:

EKG Waveforms

+ Add

EKG Comments:

Manual ECG Interpretation

+ Add

This one is for services for services who manually record ECG information

If your service does not have your cardiac monitor linked to your run form, then you will record your 12-leads in the “Manual ECG Interpretation” grid, and you must record defibrillation in the Procedures grid, and any vitals collected by your cardiac monitor (e.g. Pulse Oximetry) in the vitals grid.

Manual ECG Interpretation

+ Add Another **✓ OK** **✕ Cancel**

Date/Time of ECG:

ECG Crew Member:

ECG Type:

Method of ECG Interpretation:

Cardiac Rhythm / Electrocardiography (ECG):

Cardiac Rhythm / ECG Comments: