

~ SIREN Elite Training Highlight ~

**CARDIAC ARREST DOCUMENTATION**

Most of required documentation for a cardiac arrest call is in the “**Cardiac Arrest**” panel under the “**Patient Condition**” section. If you answer this question as “**yes**” to a cardiac arrest, the system will remind you to fill this section out with a red validation flag.

Cardiac Arrest during this incident?:

No Yes, After EMS Arrival Yes, Prior to EMS Arrival

If you **don't know the answer** to a question, you can click the “**Not Value**” square button to the right of the question. The field will be greyed out, satisfying the validation rule:

Estimated Date/Time of Cardiac Arrest: Not Known [Not Value]

If you are using a **CPR Feedback Device**, you need to document that under the question “CPR Methods”:

CPR Methods (Choose All That Apply):

Find Value...

- Device
- Compressions-External Thumper Type Device
- CPR Feedback Device**

**YOU MUST DOCUMENT ALL OF YOUR PROCEDURES AND MEDICATIONS IN THE APPROPRIATE GRIDS. It is not valid to write these in your narrative without documenting them on the form.**

To document **CPR**, go to the **Procedures** grid in the “Treatment” panel under the “Provider Action” section. Search for “**CPR**” to document one of these three options: “CPR – Cardiopulmonary Resuscitation,” “CPR with impedance threshold device,” or “Mechanical CPR Compression Device.”

Procedure:

- CPR
- CPR - Cardiopulmonary Resuscitation**
- CPR with impedance threshold device**
- Mechanical CPR Compression Device**

There are two ways to document **Defibrillation**.

<p>If you have your <b>Cardiac Monitor</b> integrated with your Run Form, a procedure will automatically be created by the monitor in the run form.</p>	<p>OR</p>	<p>If your Cardiac Monitor is <b>NOT</b> integrated with your Run Form, you must <b>manually</b> document Defibrillation in the <b>Procedures grid</b>.</p>
---	-----------	---

There are two ways to document when **No Shock is Advised**.

<p>If you have your <b>Cardiac Monitor</b> integrated with your Run Form, the "No Shock Advised" will automatically be created by the monitor in the <b>"AED/Cardiac Monitor Imports"</b> grid in the "Treatment" section.</p>	<p>OR</p>	<p>If your Cardiac Monitor is <b>NOT</b> integrated with your Run Form, you must <b>manually</b> document No Shock Advised in the <b>"Manual Cardiac Monitor Entry"</b> grid in the "Treatment" section. See below...</p>
--	-----------	---

The screenshot shows a form with the following fields and options:

- ECG Type:** A dropdown menu with "Automated External Defibrillator (AED)" selected.
- Method of ECG Interpretation:** Four buttons: "Computer Interpretation" (highlighted in blue), "Manual Interpretation", "Transmission with Remote Interpretation", and "Transmission with No Interpretation".
- Cardiac Rhythm / Electrocardiography (ECG):** A dropdown menu with "Unknown AED Non-Shockable Rhythm" selected.
- Search Bar:** A "Find Value..." search bar is located below the rhythm dropdown.

- Select "AED" for "ECG Type"
- Select "Unknown AED Non-Shockable Rhythm" for "Cardiac Rhythm"