

How to Conduct the BinaxNOW COVID-19 Antigen Self-Test

STEP 1

Obtain the test from the testing site.

STEP 2

Review all instructions associated with the test.

STEP 3

Conduct the test off-site and wait 15-30 minutes for your results.

STEP 4

Report results at healthvermont.gov/reportresults and complete the short form.

**It's important to report your test results.
Your response is confidential.**



Language Translation Services

1. Dial 1-413-216-4975
2. Enter PIN number 6187848 and “#” sign when prompted
3. Listen to automated prompt and select the language



For general inquiries, support is available from 8AM to 6PM by contacting (877) 380-3029 or support@cic-health.com

FAQs

What is a BinaxNOW™ COVID-19 Antigen Self-Test?

The BinaxNOW COVID-19 Antigen Self-Test is a lateral flow immunoassay intended for the qualitative detection of nucleocapsid protein antigen from SARS-CoV-2 from individuals with or without symptoms or other epidemiological reasons to suspect COVID-19 infection when tested twice over three days with at least 36 hours between tests. This test is authorized for non-prescription home use with self-collected direct anterior nasal (nares) swab samples from individuals aged 15 years or older, or adult-collected anterior nasal swab samples from individuals aged two years or older.

How does a BinaxNOW COVID-19 Antigen Self-Test compare to a molecular test?

Antigen tests are very specific for the virus, but are not as sensitive as molecular tests. This means that a positive result tends to be accurate, but a negative result does not rule out infection.

Why are there two tests in the kit?

Per BinaxNOW, the test is enabled for over-the-counter/self-use when people test themselves twice within three days, with 36 hours between the two tests. The two tests included allow individuals to fulfill the recommended use.

What if my test results are positive?

If you test positive, be sure to report your results via the link on other side of this card. For additional guidance, visit **healthvermont.gov/covid19positive**

To learn more about the BinaxNOW COVID-19 Antigen Self-Test, visit bit.ly/33uB6Tv